



**Symmetry Dental Equipment
Symmetry 4.0 Collection with BDS
Dual Bay Open Treatment Room
Ideal Floor Plan Layouts**


Dual Bay Open Treatment Room

The Dual Bay Open Treatment Room is designed to add an elegant feel to your practice while minimizing the amount of space required per treatment room.

Utilizing the pass-through concept from the hallway into the rear treatment consoles allows a rapid non-intrusive method to resupply the treatment room for upcoming procedures.

The Symmetry 4.0 Central Island Treatment Console reduces your startup cost by requiring only one x-ray unit per two treatment rooms.

As always we at Symmetry have designed these rooms to be standardized and fully ambidexterous making it functional for all providers. This allows for rapid growth or the sale of your practice without having to modify or purchase new equipment.



Equipment Requirements

Symmetry 4.0 4C Rear Treatment Console	2
Symmetry 4.0 SYMBDS E Delivery System and Assistant Package	2
BDS Epic Chair w FL-3000 LED Operatory Light	2
BDS DR 96 Doctor Stool	2
BDS AT 96 Assistant Stool	2
Symmetry 4.0 5A Central Island Treatment Console	1
Symmetry 4.0 3A Ancillary Treatment Console	1
Symmetry 4.0 3B Ancillary Treatment Console	1
Symmetry 4.0 3P Ancillary Treatment Console	2

Quantity

Size Requirements:

- Width: Between 19' and 22'
- Depth: Between 10' and 12'

Simple Integration with:

- Gliding monitor mount
ICW UL 500I-TS28-D3-A1
- Cabinet Mount X-Ray Console
Console has internal clearance of 52"H x 52"W x 20"D
- Flush mount nitrous panel
Plates up to 14" wide