




**Symmetry Dental Equipment
Symmetry 4.0 Collection with BDS
Single Bay Semi-Private Treatment Room
Ideal Floor Plan Layouts**

Single Bay Semi-Private Treatment Room



The Single Bay Semi-Private Treatment Room is designed to maintain patient privacy while allowing the elegant look of the Symmetry 4.0 Rear Treatment Console.

Utilizing the pass-through concept from the hallway into the rear treatment consoles allows a rapid non-intrusive method to resupply the treatment room for upcoming procedures.

The Symmetry 4.0 Ancillary Treatment Consoles maintain the elegant appearance while providing your dispensing and storage needs.

As always we at Symmetry have designed these rooms to be standardized and fully ambidexterous making it functional for all providers. This allows for rapid growth or the sale of your practice without having to modify or purchase new equipment.

Equipment Requirements

- Symmetry 4.0 4C Rear Treatment Console
- Symmetry 4.0 SYMBDS E Delivery System and Assistant Package
- BDS Epic Chair w FL-3000 LED Operatory Light
- BDS DR 96 Doctor Stool
- BDS AT 96 Assistant Stool
- Symmetry 4.0 3A Ancillary Treatment Console
- Symmetry 4.0 3B Ancillary Treatment Console
- Symmetry 4.0 3P Ancillary Treatment Console

Quantity

- 1
- 1
- 1
- 1
- 1
- 1
- 1
- 2

Size Requirements:

- Width: Between 8' 6" and 11' (may remove an ancillary console)
- Depth: Between 10' and 12'

Simple Integration with:

- Gliding monitor mount
ICW UL 500I-TS28-D3-A1
- Flush mount nitrous panel
Plates up to 14" wide