



**Symmetry Dental Equipment
Symmetry 4.0 Collection with BDS
Single Bay Surgical Suite
Ideal Floor Plan Layouts**



Single Bay Surgical Suite

The Single Bay Surgical Suite is designed to provide a very warm non-cluttered feel for your patients. Two flow patterns are available; Rear entry and front entry.

The front entry model allows the practitioners to greet the patient upon entering the room. This design allows for ample space to accommodate special needs including the use of gurneys.

The second advantage to the front entry model is the potential for internal treatment rooms with no windows. This can greatly add to the productivity of any practice.

As always we at Symmetry have designed these rooms to be standardized and fully ambidexterous making it functional for all providers. This allows for rapid growth or the sale of your practice without having to modify or purchase new equipment.



Equipment Requirements

- Symmetry 4.0 4J Rear Treatment Console
- Symmetry 4.0 SYMBDS E Delivery System and Assistant Package
- BDS Epic Chair w FL-3000 LED Operatory Light
- BDS DR 96 Doctor Stool
- BDS AT 96 Assistant Stool
- Symmetry 4.0 3A Ancillary Treatment Console
- Symmetry 4.0 3P Ancillary Treatment Console

Quantity

- 1
- 1
- 1
- 1
- 1
- 1
- 1

Size Requirements:

- Width: Between 9' 6" and 12'
- Depth: Between 10' and 12'

Simple Integration with:

- Gliding monitor mount
ICW UL 500I-TS28-D3-A1
- Flush mount nitrous panel
Plates up to 14" wide