



***Integrates with all Symmetry 4.0 Rear Treatment Consoles***



**Cabinet Sold Separately**

<p style="text-align: center;"><b>Symmetry Dental Equipment</b> <b>Symmetry 4.0 Collection</b> <b>SYMBDS A - Glide Mount Doctor Delivery Package with L/R Swing</b> <b>Specification Page</b></p>
---

**Dr. Set Up:**

- Three HP automatic
- One 3 way air / water syringe
- Wet / Dry foot control
- Silcryn tubing
- 1/4" water quick disconnect
- Easy to read pressure guage
- Pressure adjustment for each handpiece
- Non-retracting water coolant, adjustmetnt, and purge
- Master on / off toggle with indicator
- Master valve assembly with filter and regulator
- 7' umbilical tubing

**Integration Set Up:**

- Replacement base compartment door for Symmetry 4.0 series
- Cabinet mounted water bottle system

<b>Décor Finish for Door</b>	
------------------------------	--