



Integrates with all Symmetry 4.0 Rear Treatment Consoles



Cabinet Sold Separately

<p style="text-align: center;">Symmetry Dental Equipment Symmetry 4.0 Collection SYMBDS B - Glide Mount Duo Delivery Package with L/R Swing Specification Page</p>

Dr. Set Up:

- Three HP automatic
- One 3 way air / water syringe
- Wet / Dry foot control
- Silcryn tubing
- 1/4" water quick disconnect
- Easy to read pressure gauge
- Pressure adjustment for each handpiece
- Non-retracting water coolant, adjustment, and purge
- Master on / off toggle with indicator
- Master valve assembly with filter and regulator
- 7' umbilical tubing

Assistant Set Up:

- Dual autoclavable HVE
- Saliva Ejector
- One 3 way air / water syringe
- Silcryn tubing
- Vacuum canister with removable screen

Integration Set Up:

- Replacement base compartment door for Symmetry 4.0 series
- Cabinet mounted water bottle system

Décor Finish for Door	
-----------------------	--