



Integrates with all Symmetry 4.0 Rear Treatment Consoles



Cabinet Sold Separately

<p>Symmetry Dental Equipment Symmetry 4.0 Collection SYMBDS C - Glide Mount Assistant Package w L/R Swing Specification Page</p>
--

Assistant Set Up:

- Dual autoclavable HVE
- Saliva Ejector
- One 3 way air / water syringe
- Silcryn tubing
- Vacuum canister with removable screen

Integration Set Up:

- Replacement base compartment door for Symmetry 4.0 series
- Cabinet mounted water bottle system

Decor Finish for Door	
-----------------------	--