



Integrates with all Symmetry 4.0 Rear Treatment Consoles



Cabinet Sold Separately

Symmetry Dental Equipment Symmetry 4.0 Collection SYMBDS D - Arm Mount Floating Worksurface w Doctor Delivery Package Specification Page

Dr. Set Up:

- Three HP automatic
- One 3 way air / water syringe
- Wet / Dry foot control
- Silcryn tubing
- 1/4" water quick disconnect
- Easy to read pressure guage
- Pressure adjustment for each handpiece
- Non-retracting water coolant, adjustmetnt, and purge
- Master on / off toggle with indicator
- Master valve assembly with filter and regulator
- 7' umbilical tubing

Integration Set Up:

- Matching solid surface worksurface for Symmetry 4.0 series
- Cabinet mounted water bottle system

Decor Finish for Door	
Counter Finish	