



Integrates with all Symmetry 4.0 Rear Treatment Consoles



Cabinet Sold Separately

<p style="text-align: center;">Symmetry Dental Equipment Symmetry 4.0 Collection SYMBDS F - Arm Mount Floating Worksurface w Assistant Package Specification Page</p>

Assistant Set Up:

- Dual autoclavable HVE
- Saliva Ejector
- One 3 way air / water syringe
- Silcryn tubing
- Vacuum canister with removable screen

Integration Set Up:

- Matching solid surface worksurface for Symmetry 4.0 series
- Cabinet mounted water bottle system

Decor Finish for Door	
Counter Finish	